Gebrauchsinformation / Instructions for use

OdontoPlex™

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Instructions for use

OdontoPlex™

Multiplex-PCR-ELISA-System for the detection of periodontitis causing bacteria

in vitro diagnosticum

1. Introduction

The **OdontoPlex™** test system is a qualitative *in vitro* diagnostic agent for the combined molecular genetic detection of the seven periodontal pathogenic bacterial species: *Aggregatibacter actinomycetemcomitans*, *Porphyromonas gingivalis*, *Prevotella intermedia*, *Tannerella forsythia*, *Micromonas micros*, *Treponema denticola* and *Fusobacterium nucleatum*.

2. Bacteria and periodontitis

Background:

Periodontitis is a disease which can lead to the loss of teeth, if left untreated. It can also have a major effect on patients' general health. Examples of the effects of periodontitis bacteria include cardiac diseases, diseases of other internal organs and a seven fold increased risk of premature births. Unfortunately, however, the dangers of periodontitis are still considerably underestimated. And yet this disease - if diagnosed in good time - can be cured or completely prevented, if preventive measures are taken in good time. Knowledge of the factors responsible for the onset and course of the disease are of crucial importance here in selecting an individually tailored therapy. According to the latest scientific findings, genetic factors also play an important role. Newly developed test systems are able to supply the dentist with this vital information, thus paving the way at an early stage for the right prophylaxis or therapy.

Bacteria and genetic factors are responsible for the onset and course of periodontitis

Various factors are responsible for the onset, extent and course of periodontitis. That is why the course and development of periodontitis can vary from individual to individual. In order to provide a promising therapy tailored to the individual patient, it is therefore important to be able to identify the factors causing the disease in the respective patient. For example, the presence of certain types of bacteria may necessitate an antibiotic therapy. The decision whether to use an antibiotic as a therapy, and if so which antibiotic and at what dosage, depends crucially on the type and quantity of the bacteria present. In particular, the seven marker bacteria: *Aggregatibacter actinomycetemcomitans*, *Porphyromonas gingivalis*, *Tannerella forsythia*, *Prevotella intermedia*, *Micromonas micros*, *Fusobacterium nucleatum* and *Treponema denticola* are causally linked with the onset and progress of periodontitis, in that they produce and induce substances which ultimately lead to massive attachment and bone loss as a result of an inflammatory reaction.

Moreover, species of bacteria, such as *Aggregatibacter actinomycetemcomitans* and *Porphyromonas gingivalis*, have the ability to entrench themselves in the soft tissue, thereby completely evading access by instruments. Where these types of bacteria are present, a conventional, i.e. a mechanical, treatment of periodontitis is unsuccessful in the medium term, since a new bacterial colony is generally quickly established, which may possibly aggravate the patient's condition still further. Often, the only therapy option in these cases is a targeted concomitant antiobiotic therapy. It is important in this connection to select the correct and optimally effective antiobiotic in order to achieve sustained treatment success. Diagnosis before therapy using a PCR-DNA probe test is indispensable in achieving that goal.

3. Structure of the OdontoPlex™ test system

OdontoPlex™ PAF PCR reagents:

Manufacture of specific PCR amplification products from sample material, sufficient for a maximum of ten hybridisations per PCR reaction.

Suitable hybridisation modules: ***Porphyromonas gingivalis*, *Aggregatibacter actinomycetemcomitans* and *Fusobacterium nucleatum***

OdontoPlex™ PMT PCR reagents:

Manufacture of specific PCR amplification products from sample material, sufficient for a maximum of ten hybridisations per PCR reaction.

Suitable hybridisation modules: ***Tannerella forsythia*, *Prevotella intermedia*, *Micromonas micros*, *Treponema denticola* and internal control.**

OdontoPlex™ hybridisation modules:

Single wells with immobilised oligonucleotide probes for the specific detection of the relevant species (sales unit 96 wells per probe) and reagents to conduct 96 reverse hybridisations.

The following table lists the specific probes (Hyb-modules) available separately for the **OdontoPlex™** test system.

Hyb-modules for the **OdontoPlex™** test system

Symbols	Specificity	Well colour
MTS BR AA	<i>Aggregatibacter actinomycetemcomitans</i>	red
MTS BR PG	<i>Porphyromonas gingivalis</i>	blue
MTS BR FN	<i>Fusobacterium nucleatum</i>	yellow
MTS BR TF	<i>Tannerella forsythia</i>	orange
MTS BR MM	<i>Micromonas micros</i>	black
MTS BR TD	<i>Treponema denticola</i>	purple
MTS BR PI	<i>Prevotella intermedia</i>	green
MTS BR IC	Interne Kontrolle	white

4. Principle of the test

OdontoPlex™ PCR reagents contain labelled oligonucleotide primers which allow the simultaneous, specific amplification of different DNA sequences in a single PCR reaction. Specific amplification products of the different bacterial species are formed in the presence of the relevant DNA which are subsequently visualised with the different hybridisation modules.

For this purpose, the amplification products from both (**PAF** or **PMT**) PCRs are heat-denatured and added to single strand specific probes which are immobilised on the polystyrene surface of microtiter plates. A hybridisation buffer is used and a hybridisation of complementary sequences occurs which can subsequently be detected by means of the ELISA principle. For this purpose, a peroxidase (POD) conjugate is added after several stringent wash steps. This binds highly specifically to the labelling on the single strand of the PCR product bound to the oligonucleotide probe. After further wash steps, TMB substrate solution is added, which produces a blue colour after conversion by POD. This reaction is stopped by means of a stop solution. A change of colour to yellow is to be observed when this occurs. Extinction of the various wells is then measured in a photometer at a wave length of 450 nm.

Positive signals indicate the specific amplification of certain DNA sequences during the **OdontoPlex™** PCR and thus the presence of the corresponding bacterial species in the sample to be examined.

5. Reagents

5.1 Contents of PCR modules:

The reagents contained in one kit are sufficient for **96 (REF 3990)** respectively **480 (REF 3991)** determinations.

Each reagent set contains:

PRIMER PAF	Primer Mix for PAF PCR (green cap, ready for use) Contains labelled oligonucleotide primer	120 µl / 600 µl
PRIMER PMT	Primer Mix for PMT PCR (violet cap, ready for use) Contains labelled oligonucleotide primer	120 µl / 600 µl
dNTP	Nucleotide Mix (yellow cap, ready for use) Contains dATP, dCTP, dGTP, dTTP	240 µl / 1200 µl
CONTROL -	Negative control (blue cap, ready for use)	100 µl / 500 µl
MTS BR BLANK	Reagent control, breakable	16 / 48 wells

5.2 Contents of Hyb-modules:

The reagents contained in one kit are sufficient for **96 (REF 3990)** respectively **480 (REF 3991)** determinations.

Each reagent set contains:

MTS BR XX	Colour-coded wells with specific oligonucleotide probe, breakable	96 wells / 5 x 96 wells
HYBBUF	Hybridisation buffer (ready for use) Contains standard saline citrate buffer (SSC), sodium dodecyl sulphate (SDS) and N-lauroylsarcosine	60 ml / 300 ml
STRGWASH	Stringent washing solution (ready for use) Contains SSC and SDS	500 ml / 2500 ml
WASHBUF 20 X	Washing buffer (concentrated 20 times) Contains phosphate buffer, NaCl and detergent, preservatives: methylisothiazolone and oxypyron	100 ml / 500 ml
CONJ 100 X	Conjugate (concentrated 100 times , transparent cap) Preservatives: methylisothiazolone, dimethylaminoantipyrene and chloracetamide	1 ml / 5 ml
SUBS TMB	TMB (tetramethylbenzidine) substrate solution (ready for use)	100 ml / 500 ml
STOP H₃PO₄ 25%	Stop solution (25 % phosphoric acid, ready for use)	100 ml / 500 ml
	hyplex [®] Lysis buffer	60 ml / 300 ml

5.3 Additional reagents and accessories required

Reagents

- **thermostable Tth - DNA polymerase (Cat. No. 3954)**
- PCR reaction buffer
- sterile, double-distilled H₂O

Accessories

- Heating block
- PCR thermocycler
- PCR reaction vessels (recommended: 0.2 ml)
- Pipettes with sterile disposable tips: variable pipettes 0.2 - 2 µl, 2 - 20 µl and 20 - 200 µl
- Incubator (50°C ± 1°C)
- Microtiter plate photometer for measuring extinction at 450 nm and 620 nm

6. Warnings and precautions

- ☞ All reagents and materials which come into contact with potentially infectious samples must be treated with suitable disinfectant or autoclaved.
- ☞ Suitable disposable gloves must be worn during the entire test.
- ☞ The conjugate, the TMB substrate solution and the washing buffer contain preservatives. Avoid contact with the skin or mucous membranes. If this does happen, rinse affected areas with plenty of water.
- ☞ Phosphoric acid is an irritating substance. Avoid contact with the skin and mucous membranes. If this does happen, rinse affected areas with plenty of water.

7. Handling notes

The reagents of the **OdontoPlex™** PCR module should be stored upon receipt at **- 20°C** . Frequent thawing and refreezing of the reagents are to be avoided. After being used for the first time, the reagents may be stored at 2...8°C to avoid thawing and refreezing steps, but in that case they should be used within three months.

The reagents of the **OdontoPlex™** hybridisation modules should be stored at **2...8°C** .

Before the test starts, the wash buffer, TMB substrate solution, stop solution as well as the closed snap-seal bag containing the microtiter plates are to be tempered to room temperature for at least 30 minutes (18...25°C). The stringent wash solution must be preheated to 50°C. The hybridisation buffer and conjugate should always be kept cool.

The packaging bears an expiry date. Quality cannot be guaranteed after this date has passed.

In order to minimise the possibility of nonspecific product formation, it is recommended that the amplification reactions be pipetted on ice and that the DNA polymerase be the last reagent to be added. Throughout the test preparation, suitable (unpowdered) disposable gloves should be worn and pipette tips with filter insert should be used because of the danger of contamination.

8. Sample material

The sample material may consist of extracted DNA, pure bacterial cells or potentially infected material obtained from patients.

If sulcus swabs are to be tested directly with this system, 100 µl **hyplex[®] lysis buffer** per endodontic paper point is to be added to the sample transport container.

The lysis buffer is incubated for 10 minutes in a heating block (99°C). After a short centrifugation step (60 seconds at 10,000 x g), 5 µl of the supernatant is added to each of the two amplification reactions.

9. Test procedure

9.1 Preparing the amplification reaction

The reactions are prepared using the buffer supplied by the DNA polymerase manufacturer and using the following batch:

PAF – PCR:

5 µl	Sample material (e.g. 5 µl negative control)
1 µl	Nucleotide mix (yellow cap)
1 µl	PAF - Primer mix (green cap)
5 µl	10x PCR buffer
0.2 µl	Tth - DNA polymerase (1 U)
to 50 µl	double-distilled H ₂ O

and

PMT – PCR:

5 µl	Sample material (e.g. 5 µl negative control)
1 µl	Nucleotide mix (yellow cap)
1 µl	PMT - Primer mix (violet cap)
5 µl	10x PCR buffer
0.4 µl	Tth - DNA polymerase (2 U)
to 50 µl	double-distilled H ₂ O

Caution! Too much DNA or cell material may considerably reduce the effectiveness of a PCR and lead to false negative results.

9.2 Amplification reaction

The PCR thermocycler should be programmed as follows:

Cycle	Temperature [°C]	Time	Reaction
Once	94	5 min	Initial denaturation of the DNA
35 times	94	25 sec	Denaturation of the DNA
	60	25 sec	Binding of the primers
	72	20 sec + 1 sec/cycle	3'-OH extension of the primers
		45 sec	
Once	72	3 min	Final extension

After completing the PCR, the reaction mixture should be stored at 2...8°C until reverse hybridisation is conducted.

9.3 Conducting the reverse hybridisation

The eight different **OdontoPlex™** probe modules are required for conducting the reverse hybridisation. The quantities stated below refer each time to the processing of an **individual** microtiter plate well.

Examples: Five samples are to be tested for the presence of periodontal pathogenic bacteria. To do this, five probes each, **AA, PG, FN** (for the PAF-PCR) and **TF, PI, TD, MM and IC** (for the PMT-PCR) are required for the samples and 1 special probe well is required for each respective negative control. 1 reagent control well is also required. The quantities must therefore be multiplied by the factor 49 (5+5+5+5+5+5+5+5+8+1).

Before use, the stringent wash solution is preheated to 50°C and the wash buffer, TMB substrate solution, stop solution and the required microtiter plates are tempered to room temperature (18...25°C).

9.3.1 Assembling the test plates

The wells are inserted firmly into the frame. Wells not required are stored in the snap-seal bag again with the desiccant at 2...8°C.

One of the reagent controls contained in the hyplex hybridisation module must be used for each test batch to determine the reagent background signal. They are processed in the same way as the samples, after the initial addition of 50 µl hybridisation buffer. In addition, there must be a negative control for each test batch.

The PCR batch of the negative control is diluted as in section 9.3.2 and 50 µl of the solution formed is pipetted into each of the specific probe wells used and further treated in accordance with the protocol (see below).

9.3.2 Preparation of the PCR sample and hybridisation procedure

The reaction mix of the two **OdontoPlex™** PCRs carried out is denatured in the reaction vessel for 10 minutes at 95°C. It is recommended that this procedure be carried out in a thermal cycler with a heated lid.

The respective entire PCR batch (50 µl) is immediately added to **300 µl** cool (2...8°C) hybridisation solution, mixed well and **50 µl** of this is quickly transferred by pipette into each

of the corresponding wells. The microtiter plate is incubated for **30 minutes** at **50°C**. The microtiter plate should be covered in order to avoid a high evaporative loss of hybridisation solution.

9.3.3 Stringent washing

The wells are completely emptied and each are washed **briefly three times** with **200 µl** of the preheated stringent wash solution.

Care must be taken to ensure that the stringent wash solution is completely removed between the wash steps. On completion of the last wash step, tap out the plate onto a paper cloth to remove the last remnants of fluid from the wells.

9.3.4 Manufacture of the wash buffer

The wash buffer concentrate is diluted **1+19** with deionised H₂O. 0.05 ml concentrate and 0.95 ml deionised H₂O are mixed per well. The wash buffer can also be manufactured in larger quantities and stored for later use at room temperature for up to one week.

9.3.5 Washing

The wells are washed once briefly at **room temperature** with **200 µl** wash buffer.

9.3.6 Incubation with peroxidase conjugate

The conjugate solution must be prepared afresh each time.

1 µl conjugate concentrate (transparent cap) and 100 µl wash buffer per well are mixed thoroughly (dilution **1+100**) in a clean vessel. **100 µl** is added to each well by pipette and the microtiter plate is incubated for **30 minutes** at **room temperature**.

9.3.7 Washing

The wells are emptied and each is washed briefly **three times** with **200 µl** wash buffer at **room temperature**, taking care to ensure that the wash solution is completely removed between the wash steps.

9.3.8 Substrate reaction

The TMB substrate solution is ready to use. Pipette **100 µl** per well. Incubate the microtiter plate for **30 minutes** at **37°C** away from direct sunlight. The time is calculated from pipetting of the first well.

9.3.9 Stopping the reaction

To stop the reaction, pipette **100 µl** of stop solution into each well, observing the same pipetting procedure as for pipetting of the TMB substrate solution.

9.3.10 Extinction measurement

Measure the extinction of the individual wells in a microtiter plate photometer at 450 nm and a reference length of 620 – 650 nm, with zero comparison being made against air. Measurement should be carried out within 60 minutes of stopping the reaction.

10. Brief instructions on test procedure

Dilutions:

- | | | |
|---------------------------|-----------------------------------|----------------------------|
| • Washing buffer dilution | 1+19 with deion. H ₂ O | 0.05 ml + 0.95 ml per well |
| • Conjugate dilution | 1+100 with diluted washing buffer | 1 µl + 100 µl per well |

Test steps:

- | | | |
|---------------------------|---|----------------------------|
| • Denaturation | PCR reaction batch | 10 min at 95°C |
| • Sample dilution: | Add PCR reaction batch to 300 µl hybridisation solution | |
| • Hybridisation: | 50 µl per well | 30 min at 50°C ± 1° C |
| • Stringent washing: | 200 µl per well | 3 times briefly |
| • Wash step: | 200 µl per well | Once briefly |
| • Conjugate incubation: | 100 µl per well | 30 min at room temperature |
| • Wash step: | 200 µl per well | 3 times briefly |
| • Substrate incubation: | 100 µl per well | 30 min at 37°C |
| • Stop: | 100 µl per well with stop solution | |
| • Photometric measurement | 450 / 620 nm | |

11. Evaluation

The test is to be evaluated under the following conditions:

- Extinction of reagent control ≤ 0.100
- Extinction of negative control (NCE) ≤ 0.150

Evaluation of the extinction of the internal control (IC):

Extinction of the IC ≤ 0.100

Sample inhibited (or sterile) → not evaluable → new sample required;

Extinction of the IC ≤ 0.500

Sample is in principle evaluable but poor sampling, (little bacterial cell material in the sample) or patient has just received therapy;

Extinction of the IC ≥ 0.500

Sample evaluable

Evaluation scheme for direct use of patient material or DNA isolated from patients:

Evaluation criteria for PAF-PCR:

	AA	PG	FN
+	> 0.300	> 0.300	> 0.300
++	> 1.0	> 1.0	> 1.5
+++	> 2.0	> 2.0	> 2.5

Evaluation criteria for PMT-PCR:

	PI	TF	MM	TD
+	> 0.2	> 0.2	> 0.15	> 0.2
++	> 0.4	> 0.4	> 0.3	> 0.5
+++	> 0.6	> 0.6	> 0.5	> 1.0

12. Interpretation of the test results

OdontoPlex™ System is a quick, universally usable molecular biological test for detecting periodontal pathogenic bacteria in sulcus samples.

Positive results from the **OdontoPlex™** System indicate the presence of periodontal pathogenic bacteria in a sample.

The results from the **OdontoPlex™** Systems depend on the quantity and quality of the DNA in the sample.

A negative **OdontoPlex™** test result cannot rule out the presence of periodontal pathogenic bacteria in the sample material.

In addition, the effectiveness of the PCR may be reduced drastically e.g. by inhibitory substances in the sample material which may in turn result in low signals.

“Inhibited” samples should be evaluated in accordance with section 11, since sampling may have been defective.

13. Control results of the OdontoPlex™ test system

The negative control should not yield any amplification product and thus no hybridisation signal should be observable in any of the probes.

14. Troubleshooting**Constantly weak signals or no signals at all (incl. positive control)**

- Temperature of the stringent washing solution clearly above 50°C
- PCR products not (sufficiently) denatured or renatured again prior to hybridisation
- Substrate not equilibrated to room temperature
- Wrong amount of conjugate and/or substrate used
- Insufficient thermal lysis (too short; temperature not 99°C) when using the hyplex® lysis buffer

Generally high (background) signals

- The use of insufficiently heated stringent washing solution or an incubation temperature which is too low may lead to nonspecific signals. Specific, positive signals are potentiated to the same degree, thus maintaining the contrast between nonspecific and specific signals. If there is doubt as to the evaluability of the test run (see Section 11), repetition with the correct parameters is recommended.

15. Performance data

15.1 Clinical evaluation:

a) OdontoPlex™ PAF-PCR

Aggregatibacter actinomycetemcomitans

Number of PCRs	OdontoPlex™ PCR		
Culture	neg.	pos.	overall result
neg.	59	15	74
pos.	0	17	17
overall result	59	32	91

Porphyromonas gingivalis

Number of PCRs	OdontoPlex™ PCR			
Culture	bdl.*	neg.	pos.	overall result
neg.	0	46	1	47
pos.	1	2	41	44
overall result	1	48	42	91

* borderline result

Fusobacterium nucleatum

Number of PCRs	OdontoPlex™ PCR		
Culture	neg.	pos.	overall result
pos.	0	91	91
neg.	0	0	0
overall result		91	91

Evaluation of OdontoPlex™ PAF-PCR

The **sensitivities** for *Aggregatibacter actinomycetemcomitans*, *Porphyromonas gingivalis* and *Fusobacterium nucleatum* were 100 %, 95 % and 100 % respectively for the three tested detections.

The **specificities** for *Aggregatibacter actinomycetemcomitans* and *Porphyromonas gingivalis* were 79 % and 98 % respectively; no negative samples were available for *Fusobacterium nucleatum*.

b) OdontoPlex™ PMT-PCR

Treponema denticola

Number of PCRs	OdontoPlex™ PCR		
MicroIDent-PCR	neg.	pos.	overall result
neg.	2	0	2
pos.	0	8	8
overall result	2	8	10

Micromonas micros

Number of PCRs	OdontoPlex™ PCR		
Bacteriology	neg.	pos.	overall result
neg.	11	0	11
pos.	14	51	65
overall result	25	51	76

Prevotella intermedia

Number of PCRs	OdontoPlex™ PCR		
Bacteriology	neg.	pos.	overall result
neg.	33	5	38
pos.	4	34	38
overall result	37	39	76

Tannerella forsythia

Number of PCRs	OdontoPlex™ PCR		
Bacteriology	neg.	pos.	overall result
neg.	5	3	8
pos.	4	64	68
overall result	9	67	76

Evaluation of OdontoPlex™ PMT-PCR

The **sensitivities** for *Treponema denticola*, *Micromonas micros*, *Prevotella intermedia* and *Tannerella forsythensis* were 100 %, 78 %, 90 % and 94 % respectively for the four tested detections.

The **specificities** for the four pathogens were 100 %, 100 %, 87 % and 63 % (only 8 negative samples) respectively.

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