

## Safety Data Sheet

according to 2001/58/EC  
HISTO TRAY

Print date: 25.07.2007

Revision date: 25.07.2007

Page 1 of 4

---

### **1. Identification of the substance/preparation and company**

#### **1.1 Commercial product name**

HISTO TRAY B27 (Order-No. 7004, 7006, 7008, 7009), HISTO TRAY AB (Order-No. 7027, 7042, 7043), HISTO TRAY ABC (Order-No. 7013, 7022, 7035), HISTO TRAY DR (Order-No. 7029), HISTO TRAY Complement ABC (Order-No. 7003), HISTO TRAY ABC Reservoir (Order-No. 6286, 6287, 6288), HISTO TRAY AB Reservoir (Order-No. 6283, 6289), HISTO TRAY B27 Reservoir (Order-No. 6290), HISTO TRAY DR Reservoir (Order-No. 6284, 6285).

#### **Use of the substance/preparation**

For in vitro diagnostic use.

#### **1.2 Information manufacturer/supplier**

BAG Health Care GmbH

Amtsgerichtsstr. 1-5

D-35423 Lich

Telephone : +49 (0)6404 925-0

Telefax : +49 (0)6404 925-421

Contact person : Dr. med. Dietwald Dönges

Emergency telephone number : +49 (0)171-765 89 17

---

### **2. Composition/information on ingredients**

#### **Chemical characterization ( preparation )**

Microplates with predropped antisera and controls for HLA tissue typing.

#### **Additional information**

Contains no hazardous substance in concentrations requiring labelling according to EC directives currently in effect.

---

### **3. Hazards identification**

#### **Explication of special hazards for human health and environment**

No dangers to be mentioned. The Product contains material of biological origin and should there fore be handled as potentially infections. Please observe the safety information of the safety data sheet.

---

### **4. First aid measures**

#### **After inhalation**

Aerosols: Provide with fresh air. Seek medical advice.

#### **After contact with skin**

On contact with the skin, wash off immediately with soap and water.

#### **After contact with eyes**

On contact with eyes, rinse the eyes with water with the eyelids open for a sufficient length of time, then seek medical (ophthalmologic) advice immediately.

#### **After ingestion**

Rinse out mouth and then drink plenty of water. Seek medical advice.

---

## Safety Data Sheet

according to 2001/58/EC  
HISTO TRAY

Print date: 25.07.2007

Revision date: 25.07.2007

Page 2 of 4

---

### **5. Fire-fighting measures**

#### **Suitable extinguishing media**

Water. Foam. Dry powder. Carbon dioxide.

#### **Special exposure hazards arising from the substance or preparation itself, combustion products, resulting gases**

With heating up or in case of open fire poisonous gases are possible.

#### **Special protective equipment for fire fighters**

Wear self-contained breathing apparatus and protective clothing.

#### **Additional information**

Do not inhale combustible gases.

---

### **6. Accidental release measures**

#### **Personal precautions**

Ensure adequate ventilation.

#### **Environmental precautions**

Do not allow to enter drains/surface waters/groundwater. Treat recovered material as prescribed in: Disposal considerations.

#### **Methods for cleaning up/taking up**

Take up buried product. Usual absorbing materials are sufficient. Clean contaminated area with suitable disinfectants.

---

### **7. Handling and storage**

#### **7.1 Handling**

##### **Advice on safe handling**

Observe the general safety regulations for laboratories. Wear suitable protective clothing, gloves and eye/ face protection. Avoid skin and eye contact. Provide sufficient ventilation.

#### **7.2 Storage**

##### **Requirements for storage rooms and vessels**

Storage at – 20° C.

##### **Further information on storage conditions**

Consider directions for use.

---

### **8. Exposure controls/personal protection**

**Safety Data Sheet**

according to 2001/58/EC

HISTO TRAY

Print date: 25.07.2007

Revision date: 25.07.2007

Page 3 of 4

---

**8.1 Exposure limit values****Additional advice on control parameters**

not applicable

**8.2 Exposure controls****Protective and hygiene measures**

Please handle as potentially infectious. While working wear laboratory coats, eye protection and gloves. Always wash hands with water and soap after test procedure. The glove material must be impervious against the product. Due to missing tests, no recommendations regarding the selection of the glove material can be given.

---

**9. Physical and chemical properties****9.1 Appearance**

Form: fluid

Odour: odourless

**9.2 Important health, safety and environmental information**pH Value not applicable  
(at 20 °C)Water solubility: completely mixable

---

**10. Stability and reactivity****Conditions to avoid**Strong heating above room temperature.

---

**11. Toxicological information****Additional information on tests**The product has to be handled with the caution usual for chemicals.

---

**12. Ecological information****Additional information**Don't dispose to the environment.

---

**13. Disposal considerations****Advice on disposal**

All reagents, remainder and empty packing, which might contain potentially infectious material, must be inactivated before disposal (e.g. through autoclaving at 121°C). There are no standardized regulations in EEC for disposal and/ or residual substances. Chemicals, which result as residual substances, are usually hazardous waste. Their removal is regulated by appropriate laws and/ or regulations of the EEC-Member countries as well

## Safety Data Sheet

according to 2001/58/EC  
HISTO TRAY

Print date: 25.07.2007

Revision date: 25.07.2007

Page 4 of 4

---

as for the Federal Republic of Germany. Please contact the competent authority (e.g. waste disposal enterprises).

---

### **14. Transport information**

#### **14.5 Further Information**

The product is not subordinated to the transportation regulation.

---

### **15. Regulatory information**

#### **15.1 Labelling**

##### **Labelling according to EC-guidelines**

Not a hazardous product in the sense of Directives.

#### **15.2 National regulations**

Water contaminating class: - - not water contaminating  
Classification: in accordance with VwVwS appendix 4, No. 3

---

### **16. Other information**

#### **Other data**

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. The BAG Health Care GmbH shall not be held liable for any damage resulting from handling or from contact with the above product.

#### **Disclaimer**

For in vitro diagnostic use only.

---