



Material Safety Data Sheet

Product name: **Sonocheck**

updated: 09.05.05

Section 1 – Manufacture Identification

Manufacturer: PEREG GmbH
Address: Porschestraße 12
D-84478 Waldkraiburg
Germany

Manufacturer's representative: Healthmark Industries
Address: 22522 East Nine Mile road, St.Clair Shores, Mi 48080
Phone: 1-800-521-6224
Fax: 1-586-774-6473
Emergency: 1-800-521-6224

Section 2 – Product Information

Product name: SonoCheck
Product code: 30550-30

Section 3 – Hazardous Ingredients Information

Hazardous Components Of Mixture	%	OSHA PEL mg/m ³	ACGIH TLV mg/m ³	CAS Number
Bromothymol blue	<0,01	N/A	N/A	9048-46-8
Chloroform	<1	2 ppm TWA	10 ppm TWA	67-66-3

Section 4 – Physical Data

Boiling point: N/A
Solubility in water: soluble
Vapor Pressure: N/A
Appearance and odor: glass vial with green liquid

LB/GA: N/A
% volatile: N/A
Melting point: N/A

Section 5 – Fire and Explosion Hazard

Flash Point: N/A
Extinguishing media: Carbon dioxide, dry chemical, foam or water
Special Fire Fighting: N/A
Procedures: Wear appropriate protective equipment including respiratory protection as conditions warrant.

Unusual fire and Explosion hazards: Thermal decomposition may release toxic gases, such as phosgene, hydrochloric acid and chlorine
Combustion Products: Acrid smoke and irritating fumes

Flammability Limits: N/A

Material Safety Data Sheet

Sonocheck

Section 6 – Health Hazard Information

Primary Routes of

Entry:

Health Hazards: This product present no health or physical hazards under normal condition of use. **Inhalation** of decomposition products from burning may cause **eye** and **respiratory** irritation. If irritation develops following inhalation of decomposition products, remove the victim to fresh air. If symptoms persist, seek medical attention.

Section 7- First Aid measures

Eye Contact: Rinse with large amount of water . If irritation persist, seek medical advice.

Skin Contact :Rinse affected area with soap and water until no evidence of chemical remains . If irritation persists, seek medical advice.

Inhalation: N/A

Ingestion: Seek medical advise

Section 8 – Control and Protective Measures

Ventilation: not required.

Respiratory:

Protection: not required.

Skin Protection: not required.

Eye Protection: not required.

Other Protective

Measures: not required.

Section 9 – Safe Handling and Use

Spill or Release: Always wear recommended personal protective equipment. Absorb liquid with absorbing material and dispose in general waste.

Waste Disposal : General waste
Method: No special practice

This document has been prepared using data from sources considered technically reliable.

Material Safety Data Sheet

Sonocheck

It does not constitute a warranty, expressed or implied, as to the accuracy of the information contained herein. Actual conditions of use and handling are beyond the control of PEREG GmbH. Under no circumstances will PEREG GmbH be liable for the misuse or mishandling of this product.

END OF DOCUMENT

Date prepared: January 2004