




Pinpoint and fast – Precision blood group typing with real-time PCR

ERY Q Kits

BAG Diagnostics. The experts for HLA and blood group diagnostics

Your advantages




Scope of application

- Alternative to serology for rare blood groups be independent of rare blood group sera
- For donor and patient typing
- Rhesus and HNA typing during pregnancy
- For confirmatory typing of samples with weak or questionable serological results





The way it works



Coming				
SO	On! Products	Alleles	Package	REF
	Weak D-TYPE	weak D type 1, 1.1, 2, 3, 4.0/4.1, 4.2, 4.3, 5, 11, 14, 15, 17, 20, 31, 38, K409K, IVS3+1G>A	12	728401
	HPA-TYPE	HPA 1, 2, 3, 4, 5, 6, 9, 15 (a/b)	12	728402
	Partial D-TYPE	rare RHD variants (partial D)	12	728403
	HNA-TYPE	1 a/b/c, 2/2null, 3 a/a var./b, 4 a/b, 5 a/b	12	728404
	RH-TYPE	RHD/RHCE typing (C, Cw, c, D, DEL, E, e)	12	728405
	ABO-TYPE variant	ABO typing including rare variants	12	728406
	KKD / MNS-TYPE	K/k, Jk, FY, MNS	12	728407
	Rare-TYPE	Kp, Js, Di, Wr, Lu, Do, Co, Yt, Vel, Kn	12	728408

We are the experts for HLA and blood group diagnostics

With **more than 30 years of experience** BAG Diagnostics is the specialist for HLA and blood group diagnostics. As an independent, global operating company we offer **permanent ongoing development** of our products and services.

Our motivation is to **exceed the expectations of our customers** and partners. Hence, we support them with **extraordinary service** and **new ways of thinking – capable, fast and spot-on**.

BAG Diagnostics GmbH · Amtsgerichtsstr. 1-5 · 35423 Lich, Germany Tel. +496404 925-0 · Fax +496404 925-460 · info@bag-diagnostics.com · www.bag-diagnostics.com

Customer Service · Tel. +496404 925-125 · Fax +496404 925-421 · service@bag-diagnostics.com