


Liefer- und Ergebnisliste Anti-HLA-Seren 05/09


Available Anti-HLA-Sera with Test Results 05/09


0,5 ml, lyophilisiert / lyophilised




REF	Anti-HLA	LOT		€uro	Test-panel	++	f+	f-	--	%	R-Value	Bemerkungen Remarks	Menge Quantity
61099	Neg. Control	106L0384	2011-06	25,-	1013	0	0	0	1013	0			
61099	Neg. Control	106L0385	2011-11	25,-	1370	0	0	0	1370	0			
61098	Pos. Control	108L0421	2011-11	25,-	84	84	0	0	0	100			
61276	A 1	104L0263	2010-06	45,-	712	215	2	0	495	99	0,993		
61517	A 1	109L0011	2012-03	45,-	121	28	0	0	93	100	1,000		
61501	A 1+36+80	107L0653	2012-03	54,-	398	119	3	0	276	100	0,982	f+ wk A11	
61443	A 2	106L0501	2010-11	45,-	178	83	1	0	94	100	0,989		
61443	A 2	106L0502	2011-02	45,-	268	128	1	0	139	100	0,993		
61510	A 2	105L1171	2011-11	45,-	1095	524	2	0	569	100	0,996		
61371	A 2+69	105L0433	2011-07	54,-	465	216	3	0	246	100	0,987		
61328	A 2+28	105L0382	2011-06	45,-	167	84	1	0	82	100	0,988		
61428	A 3	204L0252	2011-02	45,-	2188	619	1	0	1568	100	0,999		
61428	A 3	204L0253	2011-09	45,-	2360	662	1	0	1697	100	0,999		
61518	A 3	106L0391	2012-03	45,-	1200	321	6	0	873	99	0,987		
61417	A 9 [23+24]	105L1191	2010-01	45,-	119	27	0	0	92	98	1,000	A 2403=pos.	
61417	A 9 [23+24]	105L1192	2011-11	45,-	142	31	0	0	111	99	1,000	A 2403=pos.	
61336	A 23	204L1214	2011-11	54,-	898	55	1	0	842	100	0,990		
61502	A 24	108L0161	2011-08	54,-	132	17	0	0	115	100	1,000	A 2403=pos.	
61356	A 24	304L1183	2012-03	54,-	2763	437	0	0	2326	96	1,000	A 2403=pos.	
61402	A 25	106L0303	2011-11	54,-	500	16	0	1	483	99	0,969		
61350	A 26	105L0352	2011-02	54,-	1384	92	1	0	1291	99	0,994		
REF	Anti-HLA	LOT		€uro	Test-	++	f+	f-	--	%	R-Value	Bemerkungen	Menge

					panel							Remarks	Quantity
61350	A 26	105L0353	2011-11	54,-	1558	103	1	0	1454	99	0,995		
61280	A 26+66	104L0343	2011-02	54,-	495	29	1	0	465	90	0,982		
61280	A 26+66	104L1344	2011-11	54,-	517	33	1	0	483	91	0,984		
61418	A 10 [25+26+34+66]	105L1071	2010-01	45,-	130	18	1	1	110	99	0,938	f- A 34	
61349	A 10 [25+26+34+66]	105L0453	2011-11	45,-	1198	164	4	0	1030	100	0,986	f+ A 33	
61391	A 11	106L0042	2011-09	45,-	1049	113	0	0	936	100	1,000		
61453	A 11	203L1501	2011-02	45,-	1578	155	5	0	1418	100	0,983		
61322	A 28 [68+69]	304L0875	2011-09	62,-	1949	192	1	0	1756	99	0,997		
61456	A 68+34+66	304L0732	2011-11	62,-	2282	259	41	4	1978	96	0,911	f+ A 33,25; f- A 34,66	
61429	A 29	105L1462	2010-11	62,-	687	46	5	0	636	100	0,946	f+ wk A 26	
61505	A 30	108L0371	2011-11	85,-	169	19	3	0	147	99	0,920		
61444	A 30+31	107L0961	2010-11	62,-	152	23	3	0	126	100	0,930		
61412	A 31+33 (29)	106L0413	2011-11	62,-	1214	97	34	3	1080	93	0,833	f+ A 29, A 31=s.t.wk	
61231	A 32	204L0087	2010-11	62,-	2389	170	2	2	2215	99	0,987		
61352	A 33	204L0854	2011-11	62,-	1291	43	21	0	1227	97	0,813	f+ A 10,29,31	
61430	A 33+10+28	106L0351	2010-04	45,-	303	69	19	6	209	97	0,795	f- A 28; f+ A 11	
61430	A 33+10+28	106L0352	2012-03	45,-	718	143	41	8	526	98	0,817	f- A 28, f+ A 11	
61300	A 34+25	104L0793	2011-07	62,-	544	40	1	3	500	90	0,949	s.m. A 25	
61500	B 5[51+52]	108L0151	2011-07	54,-	64	15	1	0	48	96	0,958	f+ B 53	
61500	B 5[51+52]	108L0152	2011-11	54,-	187	36	4	0	147	98	0,936	f+ B 53	
61507	B 5[51+52]+78	105L0542	2011-11	54,-	690	110	27	0	553	99	0,875	f+ A 25; B 49,53	
61383	B 5[51+52]+53	105L0572	2011-01	54,-	820	121	3	0	696	94	0,986	f+ B 49	
61383	B 5[51+52]+53	105L0574	2012-03	54,-	1203	198	3	0	1002	95	0,991	f+ B 49	
61384	B 5[51+52]+53+77	203L1223	2011-11	54,-	980	148	4	3	825	96	0,973	B 53=s.t.wk	
61511	B 5[51+52]+53+49+63	108L0291	2011-11	54,-	150	34	0	0	116	97	1,000		
61445	B 7	106L0583	2011-07	45,-	517	125	0	0	392	100	1,000		
REF	Anti-HLA	LOT		€uro	Test-panel	++	f+	f-	--	%	R-Value	Bemerkungen Remarks	Menge Quantity

61445	B 7	*	106L0584	2012-03	45,-	844	220	0	0	624	100	1,000		
61301	B 7	*	105L0043	2012-03	45,-	2271	563	8	0	1700	100	0,991	f+ B 42,55	
61446	B 7+27	*	105L1153	2012-03	45,-	210	61	4	0	145	100	0,956	f+ B 81	
61432	B 8		105L1181	2010-04	45,-	309	54	0	0	255	100	1,000		
61519	B 8	*	109L0041	2012-03	45,-	110	28	0	0	82	100	1,000		
61520	B 8	*	109L0051	2012-03	45,-	110	28	0	0	82	98	1,000		
61419	B 8+59		106L0553	2011-11	45,-	813	148	0	0	665	99	1,000		
61369	B 12 [44+45]		105L0973	2011-11	54,-	730	168	0	0	562	100	1,000		
61302	B 13		105L0032	2010-11	54,-	1300	76	2	0	1222	99	0,986		
61404	B 13	*	105L0012	2012-03	54,-	1627	99	6	0	1522	100	0,969		
61447	B 14[64+65]	*	107L0793	2012-03	54,-	566	43	0	0	523	97	1,000		
61433	B 14[64+65]		107L0062	2011-08	54,-	525	40	0	0	485	98	1,000		
61308	B 15[62+63+75+76+77]		205L0074	2011-02	54,-	1270	170	1	1	1098	98	0,993	f- B 76	
61308	B 15[62+63+75+76+77]	*	305L0075	2012-03	54,-	1692	229	2	1	1460	98	0,992	f- B 76	
61367	B 15+46+57(70)	+	204L1732	2011-11	45,-	1867	379	20	5	1463	98	0,960	f+ B 35,70; f- B 57	
61399	B 16[38+39]	*	106L0203	2012-03	54,-	1143	124	17	2	1000	96	0,921	f+ unspec.; B 38=st.wk	
61448	B 17[57+58]		106L0831	2010-11	54,-	540	52	0	0	488	96	1,000		
61413	B 17[57+58]		106L0211	2009-08	54,-	284	42	0	3	239	93	0,960	B 58 s.t.wk	
61434	B 17[57+58]+63		107L0022	2010-11	54,-	555	66	0	0	489	99	1,000		
61440	B 18		106L0452	2010-11	62,-	828	84	0	0	744	100	1,000		
61405	B 21[49+50]		105L1212	2011-11	54,-	114	11	3	0	100	100	0,873	f+ wk B 12	
61414	B 21[49+50]	+	106L0692	2011-02	54,-	592	44	2	0	546	98	0,976	f+ wk B 45	
61521	B 21[49+50]	*	107L0921	2012-03	54,-	311	23	4	0	284	95	0,917	f+ wk B 44	
61316	B 22 [54+55+56]		105L0323	2010-11	54,-	835	65	0	1	769	95	0,992	f- B 54	
61316	B 22 [54+55+56]	*	105L0325	2012-03	54,-	927	73	0	1	853	96	0,993	f- B 54	
61411	B 22[54+55+56]+42+7+67	+	106L0572	2012-03	54,-	1192	343	1	0	848	100	0,998		
REF	Anti-HLA		LOT		€uro	Test-panel	++	f+	f-	--	%	R-Value	Bemerkungen Remarks	Menge Quantity
61400	B 27		106L0232	2010-06	54,-	387	38	7	1	341	98	0,896	f+ B 47,wk B 7	

61407	B 27	105L1474	2011-11	54,-	139	13	0	0	126	100	1,000		
61332	B 27	105L0754	2011-11	54,-	1267	125	2	0	1140	99	0,991		
61333	B 27+47	105L0763	2011-11	54,-	853	120	3	1	729	97	0,981	f- B 47	
61513	B 35+53+51	108L0261	2011-11	45,-	170	56	1	2	111	99	0,961		
61468	B 35+53+5+75	105L1062	2011-11	45,-	614	198	10	0	406	100	0,964	f+ B62,78	
61436	B 37	104L0282	2010-11	62,-	2073	73	6	0	1994	100	0,960	f+ wk B44,47	
61436	B 37	204L0284	2011-11	62,-	2407	87	11	0	2309	100	0,940	f+ wk B 44,47	
61362	B 38	105L0983	2011-11	62,-	691	27	8	4	652	96	0,811	f+ B 13	
61422	B 39	93L322-2	2010-11	62,-	227	11	1	0	215	96	0,955		
61310	B 39	205L0053	2012-03	62,-	1420	65	10	0	1345	97	0,928	f+ B 67	
61415	B 39+8+14	106L0531	2009-08	54,-	173	59	1	0	113	100	0,987	f+ wk B 52	
61514	B 40[60+61]+7	108L0271	2011-11	54,-	168	61	1	0	106	99	0,987		
61449	B 40[60+61]+48(7)	204L1251	2010-11	54,-	1688	574	5	23	1086	92	0,964	f- B 7; f+ B 41	
61340	B 40[60+61]+47+13	105L0602	2011-02	54,-	1231	240	8	0	983	100	0,980		
61340	B 40[60+61]+47+13	105L0603	2012-03	54,-	1715	334	13	0	1368	100	0,977		
61359	B 41	105L0863	2011-02	85,-	346	11	3	0	332	93	0,882	s.t.wk	
61359	B 41	105L0864	2012-03	85,-	522	19	3	0	500	93	0,927	s.t.wk	
61459	B 42(7)	107L0432	2011-11	54,-	429	6	20	0	403	100	0,469	f+ B 7	
61522	B 44	108L0281	2012-03	62,-	230	51	2	0	177	99	0,976		
61425	B 44 (49)	106L0561	2010-04	62,-	156	41	4	0	111	100	0,938	f+ B 49	
61303	B 45	204L0505	2011-11	62,-	1725	30	1	0	1694	92	0,983	f+ B 82	
61319	B 45+76	94L411-2	2011-01	62,-	1980	31	5	0	1944	94	0,927	f+ B 82	
61320	B 47+44 (49)	105L0122	2011-01	62,-	1338	326	37	1	974	99	0,928	f+ B 49,52	
61364	B 49	105L0093	2010-11	62,-	1956	77	0	3	1876	91	0,980		
61450	B 49+52	107L0081	2010-11	62,-	342	20	0	1	321	93	0,974		
REF	Anti-HLA	LOT		€uro	Test-panel	++	f+	f-	--	%	R-Value	Bemerkungen Remarks	Menge Quantity
61410	B 51	106L0514	2011-11	62,-	186	20	0	0	166	97	1,000	s.t. wk	
61410	B 51	106L0515	2012-03	62,-	226	25	0	0	201	97	1,000	s.t. wk.	

61427	B 51	*	203L1843	2012-03	62,-	891	125	5	0	761	98	0,977	f+B 78, B35 homocyg.
61503	B 55		108L0171	2011-08	54,-	135	8	4	0	123	100	0,804	f+ wk B56,42
61393	B 56		106L0152	2009-09	62,-	160	5	1	0	154	100	0,910	f+ wk B55
61458	B 57		103L1171	2011-02	62,-	2771	194	16	3	2558	93	0,950	f+ wk B 58
61458	B 57		103L1172	2011-11	62,-	3037	209	17	4	2807	93	0,949	f+ wk B 58
61523	B 60	*	207L0391	2012-03	62,-	662	59	2	0	601	100	0,982	
61524	B 60	*	108L0231	2012-03	62,-	310	23	2	0	285	95	0,956	
61438	B 62		107L0202	2010-11	54,-	155	14	0	0	141	89	1,000	s.t.wk
61515	B 62		108L0111	2011-11	54,-	176	22	2	0	152	98	0,951	
61437	B 63+41+42(57)		106L0371	2010-04	62,-	766	67	2	18	679	91	0,862	f- B 57
61504	B 65+8		108L0321	2011-09	54,-	152	33	0	1	118	100	0,981	
61504	B 65+8	*	108L0322	2012-03	54,-	226	53	0	1	172	100	0,988	
61166	B 70[71+72]+15+21+56		94L402-4	2010-11	62,-	1802	475	108	7	1212	96	0,855	f+ B 35,52,57
61312	B 70+62+35+45+41+50(60)		94L403-4	2011-11	62,-	2907	1167	8	19	1713	95	0,981	f- B 60,62; f+ B 44,49,75
61466	B 73	*	205L1502	2012-03	62,-	1665	7	2	0	1656	93	0,881	
61394	Bw 4+B 45		105L0173	2011-04	54,-	911	538	1	36	336	92	0,918	f- B51,17,37,38
61394	Bw 4+B 45	*	105L0175	2012-03	54,-	973	566	1	38	368	93	0,919	f- B51,17,37,38
61516	Bw 6		108L0431	2011-11	54,-	63	47	0	8	8	100	0,654	f- B 35,18,70
61516	Bw 6	*	108L0432	2012-03	54,-	170	131	2	12	25	99	0,746	f- B 35,18,70
61506	Bw 6	*	108L0241	2011-11	54,-	237	194	5	5	33	98	0,843	f- B 35; f+ B 13

REF	Anti-HLA	LOT		€uro	Test-panel	++	f+	f-	--	%	R-Value	Bemerkungen Remarks	Menge Quantity	
61408	Cw 1	106L0312	2010-04	54,-	669	51	2	0	616	97	0,979			
61408	Cw 1	*	106L0314	2012-03	54,-	1149	85	2	0	1062	98	0,988		
61409	Cw 2	106L0293	2011-07	54,-	901	82	0	0	819	95	1,000			

61409	Cw 2	*	106L0294	2012-03	54,-	1151	106	0	0	1045	96	1,000		
61293	Cw 3		94L385-3	2009-08	54,-	591	158	20	2	411	98	0,911	f+ unspec.	
61470	Cw 3	*	107L0702	2012-03	54,-	520	120	0	0	400	99	1,000		
61452	Cw 4		107L0881	2011-01	54,-	143	20	0	0	123	100	1,000		
61365	Cw 5		204L1232	2010-04	54,-	1277	168	0	2	1107	92	0,993	s.t.weak	
61471	Cw 5		107L0442	2011-11	54,-	342	33	1	0	308	99	0,984		
61342	Cw 6		94L428-4	2011-11	85,-	1294	247	35	0	1012	98	0,920	f+ wk Cw 4	
61398	Cw 7		106L0182	2011-09	54,-	306	156	3	9	138	90	0,922		
61398	Cw 7	*	106L0183	2012-03	54,-	380	200	6	10	164	91	0,915		
61509	Cw 8+5	*	104L1422	2011-11	62,-	129	19	0	2	108	93	0,943		
5931	Pos.B-Cell-Control		93L908-1	2011-03	30,-	20	20	0	0	0	100			
5930	Pos.T-Cell-Control		106L0811	2011-03	30,-	80	80	0	0	0	100			

s.t. = manchmal/*some times*; wk = schwach/*weak*; * = Neu im Lieferprogramm/*New specificities and antisera*; + = geringe Stückzahl verfügbar / *small amount available*

Ort/Datum
Date

Lieferanschrift/Unterschrift
Address/Signature