

CERTIFICATE

BAG-BlackCheck C

It is herewith certified that the chemical indicator BAG-BlackCheck C (REF 7514), from BAG Health Care GmbH fulfils the requirements of a **Class 4** indicator according to the International Norm **ISO 11140-1**.

The certification is limited to the use in steam sterilization processes at a sterilization temperature of 121 °C and 134 °C respectively.

The product presents no health or physical hazards under normal conditions of use.

BAG Health Care GmbH

August 2007



Dr. Lenard Mueller
Product Manager
Hygiene Monitoring