

# CERTIFICATE

## BAG-ChemoStrip 134

It is herewith certified that the chemical indicator BAG-ChemoStrip 134 (REF 7502), from BAG Health Care GmbH fulfils the requirements of an integrating **Class 6** indicator according to the International Norm **ISO 11140-1**.

The certification is limited to the use in steam sterilization processes at a sterilization temperature of 134 °C.

The chemical indicator BAG-ChemoStrip 134 is free of chlorine, lead, cadmium and thallium. The product presents no health or physical hazards under normal conditions of use.

BAG Health Care GmbH

August 2007



Dr. Lenard Mueller  
Product Manager  
Hygiene Monitoring