

# CERTIFICATE

## BAG-HeliPac

It is herewith certified that the chemical indicators BAG-HeliPac 121° C (REF 7571), BAG-HeliPac 134° C (REF 7572) und BAG-HeliPac 18 (REF 7573) from BAG Health Care GmbH fulfils the requirements of an emulating **Class 6** indicator according to the International Norm **ISO 11140-1**.

The certification is limited to the use in steam sterilization processes at a sterilization temperature of 121° C and 134° C respectively.

The chemical indicators BAG-HeliPac is free of chlorine, lead, cadmium and thallium. The product presents no health or physical hazards under normal conditions of use.

BAG Health Care GmbH

August 2007



Dr. Lenard Mueller  
Product Manager  
Hygiene Monitoring