

CERTIFICATE

BAG-Integraph

It is herewith certified that the chemical indicator BAG-Integraph (REF 7510), from BAG Health Care GmbH fulfils the requirements of an integrating **Class 5** indicator according to the International Norm **ISO 11140-1**.

The certification is limited to the use in steam sterilization processes at a sterilization temperature of 121 °C and 134 °C respectively.

The chemical indicator BAG-Integraph ist free of chlorine, lead, cadmiun and thallium. The product presents no health or physical hazards under normal conditions of use.

BAG Health Care GmbH

August 2007



Dr. Lenard Mueller
Product Manager
Hygiene Monitoring